
(DE)	Gebrauchsanweisung .....	3
(GB)	Instructions for use .....	6
(FR)	Mode d'emploi .....	9
(ES)	Instrucciones de uso .....	12
(IT)	Istruzioni per l'uso .....	15

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Stand: 04.2009

**Graphite-free-phosphate-bonded, precision investment material for investing all types of chrome cobalt molybdenum and in precious metal - for use with silicone and gel duplicating materials**

### Storage

Keep the powder in a cool, dry place.

Do not store the mixing liquid below +5 °C as it is sensitive to low temperatures.

### Working temperature

Approx. 22 °C (room temperature)

### Mixing ratio:

#### Silicone-Duplicating

100 g Powder – 18 ml Liquid

160 g Powder – 29 ml Liquid

#### Gel-Duplicating

100 g Powder – 16 to 18 ml Liquid

160 g Powder – 25 to 29 ml Liquid

If the amount of liquid is reduced, the concentration must be adapted accordingly, e.g. to 50 % for 160 g/25 ml

Liquid concentration for 1000 ml Liquid		
Liquid concentration (examples)		
Concentration	Concentrate (ml)	dist. water (ml)
50 %	500	500
65 %	650	350
75 %	750	250

Concentrations							
		Sachet 1 x 160 g Fabricating 1 Model		Sachet 2 x 160 g (320 g) Fabricating 2 Models		Sachet 3 x 160 g (480 g) Fabricating 3 Models Investing 1 Ring	
Mixing ratio Powder: Liquid	Concen- tration %	Liquid (ml)	dist. water (ml)	Liquid (ml)	dist. water (ml)	Liquid (ml)	dist. water (ml)
Silicone-Duplicating							
160 g – 29 ml	75	22	7	43	15	65	22
	65	19	10	37	21	56	31
Gel-Duplicating *							
160 g – 29 ml	75	22	7	43	15	65	22
	65	19	10	37	21	56	31
160 g – 27 ml	65	17	10	35	19	53	28
	50	13	14	27	27	40	41
160 g – 25 ml	50	12	13	25	25	37	38

If the amount of liquid is reduced, the concentration must be adapted accordingly, e.g. to 50 % for 160 g/25 ml

\* The model must only be duplicated with gel if it has been fully invested in investment

Expansion control (%)				
	Silicone-Duplicating		Gel-Duplicating *	
	Model	Investing	Model	Investing
Cobalt-chrome-Molybdenum (CoCrMo)	75	75	50–75	50–75
Gold casting alloy	65	65	65	50–65
*also refer to the mixing ratio				

### **Stirring**

Use a spatula to stir the investment material by hand until the powder has been wetted thoroughly.

### **Mixing in vacuum**

60 seconds.

Place in vacuum for 15 seconds but without switching the paddles on.

Mix for 60 seconds. Check the vacuum/paddles frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the casting.

### **Working time**

Approx. 4 minutes at approx. 22 °C (room temperature)

The working time depends on the room temperature. Warmth curtails the working time.

## **Fabricating the model in duplicating silicone**

### **Fabricating investment model**

Vibrator speed 6000 min<sup>-1</sup>

Before casting the silicone duplicating mould, spray it with Hera SWE 2000 silicone wetting agent, wait 2–3 minutes and blow it dry. The excellent flow properties facilitate casting non-porous models. Set the vibrator to high frequency and low intensity. Once the duplicating mould has been covered, switch the vibrator off immediately. Vibrating too intensely causes the material to segregate, impairs the precision of fit and results in rough surfaces (bubbles) as well as sedimentation.

### **Drying time Investment model**

30 minutes. Carefully remove the model after 30 minutes. We recommend applying adhesive to the model or heating the working model briefly to approx. 50 °C to ensure that the wax pattern adheres properly.

## **Fabricating the model in duplicating gel**

### **Fabricating investment model**

Vibrator speed 6000 min<sup>-1</sup>

Carefully remove the plaster model from the gel duplicating mould. Ensure that the mould is not distorted. The excellent flow properties facilitate casting non-porous models. Set the vibrator to high frequency and low intensity. Once the duplicating mould has been covered, switch the vibrator off immediately. Vibrating too intensely causes the material to segregate, impairs the precision of fit and results in rough surfaces (bubbles) as well as sedimentation.

### **Drying time Investment model**

30 minutes / Carefully remove the model after 30 minutes.

### **Time taken for model to dry**

60 minutes at 220 °C

### **Hardening the model**

The model should be a light colour before immersion in the cold hardening liquid for 5–8 seconds. After immersion, return the model to the drying chamber for 2 minutes.

### **Full flask investing**

Vibrator speed 6000 min<sup>-1</sup>

Set the vibrator to high frequency and low intensity. Switch off the vibrator as soon as the pattern is covered with investment to avoid possible segregation of the investment (sedimentation). The model must only be duplicated with gel if it has been fully invested in investment.

**The model must only be duplicated with gel if it has been fully invested in investment.**

Suitable ring systems
a) Heraeus Kulzer metal denture base system with paper sleeve
b) Other alternatives can also be used for casting metal denture bases

#### Placing the rings in a preheated furnace

Allow the investment mould to set for at least 30 minutes before placing it in the pre-heating furnace. Roughen or grind the top and bottom surfaces of the investment mould before placing it in the furnace.

#### Preheating

Preheating process for CoCrMo frameworks			
Program level	Temperature (°C)	Heat-rate (°C/min.)	Hold time (min.)
1	e.g. 23 °C room temp. – 270 °C	8	60
2	270 °C – 580 °C	8	60
3	580 °C – 1030 °C	8	40

Preheating process for precious metal alloy frameworks			
Program level	Temperature (°C)	Heat-rate (°C/min.)	Hold time (min.)
1	e.g. 23 °C room temp. – 270 °C	8	60
2	270 °C – 580 °C	8	60
3	580 °C – 700 °C	8	40

#### ⚠ Safety Instructions

- Do not open the furnace during preheating as this would cause aggressive combustion.
- The preheat furnace should only be filled to half of its maximum capacity.
- Investments contain quartz. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer). Recommendation: Wear a protective face mask, type FFP 2 EN 149:2001. Prevent the release of dust when opening the bag and filling the mixing bowl.
- Rinse the empty bag with water before crumpling it up.
- Always moisten any dust before removing it from the workplace.
- After casting, place the completely cooled mould in water until it is fully moistened to prevent dust when devesting.
- Use an extractor with a fine dust filter when sandblasting.

We always attempt to include the latest results for our investment materials (development and application) in the instructions. Excellent precision of fit and surfaces are achieved with Heraeus equipment and materials. The results may vary if other materials (e.g. casting ring liners, waxes, equipment etc.) or different techniques are used.

Should you have any further questions, please call our product specialists.

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**Contact in Australia****Heraeus Kulzer Australia Pty. Ltd.**

Locked Bag 750  
Roseville NSW 2069  
Phone +61 29.417-8411  
Fax +61 29.417-5093  
sales@heraeus.com.au  
www.heraeus-dental.com

**Contact in France****Heraeus, Division Dentaire**

12, Avenue du Québec  
Villebon - B.P.630  
91945 Courtabœuf Cédex  
Phone +33 169.18.48.85  
Fax +33 169.28.78.22

dentaire@heraeus.com  
www.heraeus-dental.com

**Contact in Scandinavia****Heraeus Kulzer Nordic AB**

Box 437  
SE-191 24 Sollentuna  
Hammarbacken 4B  
Phone +46 8585.777.55  
Fax +46 8623.14.13  
nordic.dental@heraeus.com  
www.heraeus.se

**Contact in Mexico****Heraeus Kulzer Mexico S.A. de C.V.**

Homero 527  
301 y 302 Co. Pol.  
11560 Mexico  
Phone +52 55.5531-5549  
Fax +52 55.5525-1651  
contacto@dheraeus-kulzer.com.mx  
www.heraeus-kulzer.com.mx

**Contact in Benelux**

Heraeus Kulzer Benelux B.V.  
Postbus 986  
NL-2003 RZ Haarlem  
Phone +31 23.543.42-50  
Fax +31 23.543.42-55  
info-benelux@heraeus.com  
www.heraeus-dental.com

**Contact in Italy****Heraeus Kulzer S.r.l.**

Via Console Flaminio 5/7  
20134 Milano  
Phone +39 02210.09.41  
Fax +39 02210.09.42-83  
heraeus.hki@heraeus.com  
www.heraeuskulzer.it

**Contact in Spain****Heraeus S.A.**

Forjadores, 16  
Prado del Espino  
28660 Boadilla del Monte  
Madrid  
Phone +34 91358.03-75  
Fax +34 91358.03-68  
dental@heraeus.es  
www.heraeus-dental.com

**Contact in Austria****Heraeus Kulzer Austria GmbH**

Nordbahnstr. 36/2/4/ Top 4.5  
1020 Wien  
Phone +43 1.408.09.41  
Fax +43 1.408.09.41-70  
officehkat@heraeus.com  
www.heraeus-kulzer.at

**Contact in Switzerland****Heraeus Kulzer Switzerland AG**

Ringstrasse 15A  
8600 Dübendorf  
Phone +41 43.333.72-50  
Fax +41 43.333.72-51  
officehkch@heraeus.com  
www.heraeus-kulzer.ch

**Contact in the United Kingdom****Heraeus Kulzer Ltd.**

Heraeus House, Albert Road  
Northbrook Street, Newbury  
Berkshire, RG14 1DL  
Phone +44 1635.30-500  
Fax +44 1635.30-606  
sales@kulzer.uk  
www.heraeus-dental.com

Heraeus Kulzer GmbH  
Grüner Weg 11  
63450 Hanau  
Phone 0800 43 72-522  
info.lab@heraeus.com  
www.heraeus-dental.com

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